

## Method Information

GC-ECD  
MethodECD.M  
OC/OP/ON Scan by Steven C. Moser  
Environmental/Chemical Laboratory Scientist

GC-uECD/EPC

ODAFF Pesticide Laboratory Residue Section

## Method Change History

Operator	Date	Change Information
Keith Keesee	4/19/2012 7:50:04 AM	This method was created at 4/19/2012 7:50:04 AM and based on method C:\HPCHEM\1\METHODS\SCM-135.M
Keith Keesee	4/19/2012 7:50:05 AM	
Keith Keesee	4/19/2012 7:50:21 AM	
Steven C. Moser	4/24/2012 11:32:34 AM	This method was created at 4/24/2012 11:32:34 AM and based on method C:\HPCHEM\1\METHODS\SCM-60.M
Steven C. Moser	4/24/2012 11:32:41 AM	
Steven C. Moser	4/24/2012 11:34:20 AM	
Keith Keesee	5/2/2012 9:18:09 AM	
Keith Keesee	5/2/2012 9:44:21 AM	Adj. Integration Events
Keith Keesee	5/2/2012 9:49:15 AM	Adj. Peak Labels
Keith Keesee	5/2/2012 9:52:55 AM	
Keith Keesee	5/2/2012 11:25:47 AM	This method was created at 5/2/2012 11:25:47 AM and based on method C:\HPCHEM\1\METHODS\ECD.M
Keith Keesee	5/2/2012 11:25:48 AM	
Keith Keesee	5/2/2012 11:39:04 AM	
Keith Keesee	5/2/2012 5:31:22 PM	Adj. Integration Events
Keith Keesee	5/3/2012 3:22:23 PM	Adj. RTs
Keith Keesee	5/3/2012 3:23:49 PM	Adj. RTs
Keith Keesee	5/3/2012 5:36:21 PM	This method was created at 5/3/2012 5:36:21 PM and based on method C:\HPCHEM\1\METHODS\ECD2.M
Keith Keesee	5/3/2012 5:36:23 PM	
Keith Keesee	5/23/2012 12:20:58 PM	Adj. RT
Keith Keesee	6/12/2012 11:37:40 AM	Update EPA 508 RTs
Keith Keesee	6/12/2012 11:50:40 AM	
Keith Keesee	6/12/2012 11:52:32 AM	
Keith Keesee	6/12/2012 3:49:14 PM	
Keith Keesee	6/12/2012 3:54:03 PM	
Keith Keesee	6/12/2012 3:56:47 PM	
Keith Keesee	6/12/2012 4:10:51 PM	
Keith Keesee	6/13/2012 11:18:46 AM	
Keith Keesee	6/13/2012 11:21:14 AM	
Keith Keesee	6/15/2012 9:40:37 AM	
Keith Keesee	6/15/2012 9:47:23 AM	
Keith Keesee	6/15/2012 9:58:14 AM	adjusted Intergration Prams
Keith Keesee	6/15/2012 10:05:02 AM	
Keith Keesee	6/15/2012 10:21:00 AM	updated Calbration Table

## Method Change History

Operator	Date	Change Information
Keith Keesee	7/12/2012 9:16:40 AM	Adj. RT of Permethrin
Keith Keesee	8/2/2012 11:48:11 AM	
Keith Keesee	8/14/2012 9:52:15 AM	
Keith Keesee	8/23/2012 1:03:44 PM	
Keith Keesee	10/24/2012 4:27:19 PM	
elyon	11/1/2012 9:24:52 AM	
Jennifer Bledsoe	5/10/2013 2:33:53 PM	Update RT
Elena Lyon	5/16/2013 1:13:17 PM	
Jennifer Bledsoe	5/30/2013 3:57:24 PM	
Jennifer Bledsoe	6/14/2013 4:08:26 PM	Installed new column Serial #1086027
Jennifer Bledsoe	6/14/2013 4:09:48 PM	
Jennifer Bledsoe	6/17/2013 9:39:12 AM	
Jennifer Bledsoe	6/17/2013 10:01:28 AM	
Jennifer Bledsoe	6/17/2013 4:08:43 PM	
Jennifer Bledsoe	6/18/2013 3:21:31 PM	Updated RT for InCk and EPA Mix after baking column
Jennifer Bledsoe	6/19/2013 9:33:58 AM	
Jennifer Bledsoe	6/27/2013 8:19:37 AM	Updated RTs
Jennifer Bledsoe	6/28/2013 8:10:38 AM	
Jennifer Bledsoe	6/28/2013 9:19:12 AM	Updated RT
Jennifer Bledsoe	6/28/2013 11:47:41 AM	
Jennifer Bledsoe	7/23/2013 3:28:28 PM	This method was created at 7/23/2013 3:28:28 PM and based on method C:\HPCHEM\1\DATA\070213E1\CAPTAN.D\RUN.M
Jennifer Bledsoe	7/23/2013 3:28:30 PM	
Jennifer Bledsoe	8/8/2013 8:27:55 AM	
Jennifer Bledsoe	8/8/2013 8:38:09 AM	
Jennifer Bledsoe	8/8/2013 8:42:44 AM	
Jennifer Bledsoe	8/8/2013 8:47:06 AM	
Jennifer Bledsoe	8/8/2013 8:51:54 AM	
Jennifer Bledsoe	8/8/2013 8:58:04 AM	
Jennifer Bledsoe	8/8/2013 1:47:21 PM	Updated EPA mix cal curve
Laura Esquivel	10/29/2013 3:20:15 PM	
Laura Esquivel	10/30/2013 3:06:16 PM	
Kyle Baker	11/1/2013 12:46:15 PM	
Kyle Baker	11/8/2013 10:09:14 AM	
Kyle Baker	11/15/2013 9:38:39 AM	
Kyle Baker	11/15/2013 10:19:19 AM	
Kyle Baker	11/18/2013 11:48:02 AM	
Kyle Baker	11/18/2013 11:53:32 AM	
Kyle Baker	12/11/2013 11:12:04 AM	
Kyle Baker	12/11/2013 11:14:09 AM	
Kyle Baker	3/18/2014 12:16:01 PM	This method was created at 3/18/2014 12:16:01 PM and based on method C:\HPCHEM\1\DATA\010614E2\505.D\RUN.M
Kyle Baker	3/18/2014 12:16:16 PM	Had to fix method after a power bump. SM
Kyle Baker	3/18/2014 12:17:22 PM	
Kyle Baker	5/5/2014 9:11:52 AM	
Laura Esquivel	5/9/2014 8:43:12 AM	
Kyle Baker	6/23/2014 10:52:37 AM	set split to 25:1 instead of 5:1

Method Change History

Operator	Date	Change Information
Kyle Baker	6/24/2014 3:43:33 PM	Changed split ratio back to 5:1
Jennifer Busey	7/1/2014 3:22:19 PM	
Jennifer Busey	7/2/2014 8:43:59 AM	Updated RT
Jennifer Busey	8/4/2014 1:57:07 PM	
Jennifer Busey	8/14/2014 10:03:42 AM	RT Updated
Jennifer Busey	8/26/2014 8:47:17 AM	
Jennifer Busey	8/26/2014 8:48:25 AM	integrated Diazinon
Yvette Turner	9/24/2014 12:38:32 PM	
Yvette Turner	9/24/2014 12:41:14 PM	
Yvette Turner	10/7/2014 9:16:53 AM	peak integration
Yvette Turner	10/7/2014 9:46:56 AM	
Yvette Turner	10/7/2014 11:15:07 AM	
Yvette Turner	10/7/2014 11:32:22 AM	
Yvette Turner	11/4/2014 8:00:57 AM	update retention times
Yvette Turner	11/4/2014 8:44:00 AM	
Yvette Turner	11/4/2014 9:13:46 AM	
Yvette Turner	2/12/2015 9:26:00 AM	
Yvette Turner	2/12/2015 11:05:31 AM	updated retention times
Yvette Turner	2/12/2015 12:09:54 PM	
Yvette Turner	2/17/2015 8:53:19 AM	Updated retention times
Yvette Turner	3/9/2015 10:44:51 AM	This method was created at 3/9/2015 10:44:51 AM and based on method C:\HPCHEM\1\DATA\022615E1\345.D\RUN.M
Yvette Turner	3/9/2015 10:44:52 AM	

Run Time Checklist

Pre-Run Cmd/Macro: off  
 Data Acquisition: on  
 Standard Data Analysis: on  
 Customized Data Analysis: off  
 Save GLP Data: on  
 Post-Run Cmd/Macro: off  
 Save Method with Data: skipped - no ACQ running

Injection Source and Location

Injection Source: GC Injector  
 Injection Location: Front

=====  
 6890 GC METHOD  
 =====

## OVEN

Initial temp: 180 'C (On)                   Maximum temp: 350 'C  
 Initial time: 1.50 min                   Equilibration time: 0.50 min

## Ramps:

#	Rate	Final temp	Final time
1	30.00	220	2.17
2	30.00	280	0.00
3	40.00	300	5.00
4	0.0(Off)		

Post temp: 180 'C  
 Post time: 0.00 min  
 Run time: 12.50 min

## FRONT INLET (SPLIT/SPLITLESS)

Mode: Pulsed Split  
 Initial temp: 230 'C (On)  
 Pressure: 9.57 psi (On)  
 Split ratio: 5:1  
 Pulse pressure: 20.0 psi  
 Pulse time: 0.50 min  
 Split flow: 12.5 mL/min  
 Total flow: 21.4 mL/min  
 Gas saver: On  
 Saver flow: 15.0 mL/min  
 Saver time: 60.00 min  
 Gas type: Hydrogen

## BACK INLET (SPLIT/SPLITLESS)

Mode: Split  
 Initial temp: 210 'C (On)  
 Pressure: 5.94 psi (On)  
 Split ratio: 5.002:1  
 Split flow: 7.1 mL/min  
 Total flow: 15.9 mL/min  
 Gas saver: Off  
 Gas type: Hydrogen

## COLUMN 1

Capillary Column  
 Model Number: Restek 13324  
 Rxi-MS  
 Max temperature: 350 'C  
 Nominal length: 30.0 m  
 Nominal diameter: 320.00 um  
 Nominal film thickness: 0.25 um  
 Mode: ramped flow  
 Initial flow: 2.5 mL/min  
 Initial time: 4.00 min

#	Rate	Final flow	Final time
1	4.40	18.0	0.00
2	9.40	28.0	5.50
3	0.0(Off)		

Post flow: 0.0 mL/min  
 Nominal init pressure: 9.59 psi  
 Average velocity: 59 cm/sec  
 Inlet: Front Inlet  
 Outlet: Front Detector  
 Outlet pressure: ambient

## COLUMN 2

Capillary Column  
 Model Number: Restek 13324  
 Rxi-1MS, 30m x 0.32mm x 0.25um  
 Max temperature: 350 'C  
 Nominal length: 30.0 m  
 Nominal diameter: 320.00 um  
 Nominal film thickness: 0.25 um  
 Mode: constant pressure  
 Pressure: 5.94 psi  
 Nominal initial flow: 1.4 mL/min  
 Average velocity: 37 cm/sec  
 Inlet: Back Inlet  
 Outlet: Back Detector  
 Outlet pressure: ambient

## FRONT DETECTOR (μECD)

Temperature: 325 'C (On)  
 Mode: Constant makeup flow  
 Makeup flow: 135.0 mL/min (On)  
 Makeup Gas Type: Nitrogen  
 Electrometer: On

## BACK DETECTOR (NPD)

Temperature: 325 'C (On)  
 Hydrogen flow: 2.0 mL/min (Off)  
 Air flow: 60.0 mL/min (Off)  
 Mode: Constant column+makeup flow  
 Combined flow: 15.0 mL/min  
 Makeup flow: On  
 Makeup Gas Type: Nitrogen  
 Adjust offset: 30.00  
 Electrometer: On  
 Bead: Off  
 Equilibration time: 5.00

SIGNAL 1

Data rate: 50 Hz  
Type: front detector  
Save Data: On  
Zero: 0.0 (Off)  
Range: 0  
Fast Peaks: Off  
Attenuation: 0

SIGNAL 2

Data rate: 20 Hz  
Type: back detector  
Save Data: Off  
Zero: 0.0 (Off)  
Range: 0  
Fast Peaks: Off  
Attenuation: 0

COLUMN COMP 1

Derive from front detector

COLUMN COMP 2

Derive from back detector

POST RUN

Post Time: 0.00 min

TIME TABLE

Time	Specifier	Parameter & Setpoint
------	-----------	----------------------

GC Injector

Front Injector:

Sample Washes	2
Sample Pumps	3
Injection Volume	1.0 microliters
Syringe Size	10.0 microliters
PostInj Solvent A Washes	3
PostInj Solvent B Washes	1
Viscosity Delay	0 seconds
Plunger Speed	Fast
PreInjection Dwell	0.00 minutes
PostInjection Dwell	0.00 minutes
Sampling Depth	0.0 mm

Back Injector:

No parameters specified

=====  
 Integration Events  
 =====

Results will be produced with the enhanced integrator.

-----  
 Default Integration Event Table "Event"  
 -----

Event	Value	Time
Initial Slope Sensitivity	1.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.700	Initial
Initial Shoulders	OFF	Initial

-----  
 Detector Default Integration Event Table "Event\_TCD"  
 -----

Event	Value	Time
Initial Slope Sensitivity	100.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.000	Initial
Initial Shoulders	OFF	Initial

-----  
 Detector Default Integration Event Table "Event\_ADC"  
 -----

Event	Value	Time
Initial Slope Sensitivity	20.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.000	Initial
Initial Shoulders	OFF	Initial

-----  
 Detector Default Integration Event Table "Event\_FID"  
 -----

Event	Value	Time
Initial Slope Sensitivity	50.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.000	Initial
Initial Shoulders	OFF	Initial

-----  
 Detector Default Integration Event Table "Event\_FPD"  
 -----

Event	Value	Time
Initial Slope Sensitivity	50.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.000	Initial
Initial Shoulders	OFF	Initial

-----  
Detector Default Integration Event Table "Event\_uECD"  
-----

Event	Value	Time
Initial Slope Sensitivity	146.460	Initial
Initial Peak Width	0.030	Initial
Initial Area Reject	12.000	Initial
Initial Height Reject	6.000	Initial
Initial Shoulders	OFF	Initial

-----  
Detector Default Integration Event Table "Event\_ECD"  
-----

Event	Value	Time
Initial Slope Sensitivity	100.000	Initial
Initial Peak Width	0.080	Initial
Initial Area Reject	1.000	Initial
Initial Height Reject	1.000	Initial
Initial Shoulders	OFF	Initial

-----  
Detector Default Integration Event Table "Event\_NPD"  
-----

Event	Value	Time
Initial Slope Sensitivity	500.000	Initial
Initial Peak Width	0.040	Initial
Initial Area Reject	0.010	Initial
Initial Height Reject	0.100	Initial
Initial Shoulders	OFF	Initial

-----  
Signal Specific Integration Event Table "Event\_ECD1A"  
-----

Event	Value	Time
Initial Slope Sensitivity	2000.000	Initial
Initial Peak Width	0.005	Initial
Initial Area Reject	100.000	Initial
Initial Height Reject	50.000	Initial
Initial Shoulders	DROP	Initial
Integration	OFF	0.000
Integration	ON	1.000
Integration	OFF	9.000

Apply Manual Integration Events: No

Advanced Baseline : Yes

=====  
Calibration Table  
=====

Calib. Data Modified : 2/17/2015 8:50:10 AM

Calculate : External Standard  
Based on : Peak HeightRel. Reference Window : 1.000 %  
Abs. Reference Window : 0.000 min

Rel. Non-ref. Window : 1.000 %  
 Abs. Non-ref. Window : 0.000 min  
 Do not use Multiplier & Dilution Factor with ISTDs  
 Uncalibrated Peaks : compound name not specified  
 Partial Calibration : Yes, identified peaks are recalibrated  
 Correct All Ret. Times: Yes, even for non-identified peaks

Curve Type : Linear (some peaks differ, see below)  
 Origin : Ignored (some peaks differ, see below)  
 Weight : Equal (some peaks differ, see below)

Recalibration Settings:  
 Average Response : Average all calibrations  
 Average Retention Time: No Update

Calibration Report Options :  
 Printout of recalibrations within a sequence:  
     Calibration Table after Recalibration  
     Normal Report after Recalibration  
 If the sequence is done with bracketing:  
     Results of first cycle (ending previous bracket)

Signal 1: ECD1 A,

RetTime [min]	Lvl Sig	Amount [ng/ul]	Height	Amt/Height	Ref Grp Name
1.027	1	4.42000e-1	1424.56201	3.10271e-4	Clopyralid m.e.
1.295	1	2.12000e-1	439.90289	4.81925e-4	Dicamba m.e.
1.366	1	1.00000	1945.69043	5.13956e-4	Methomyl
1.710	1	4.48000e-1	313.02985	1.43117e-3	2,4-D
1.812	1	1.00000e-1	5868.43701	1.70403e-5	TCMX
2.058	1	1.00000e-2	4696.91211	2.12906e-6	alpha-BHC
		5.00000e-2	6.63500e4	7.53580e-7	
		1.00000e-1	1.86725e5	5.35547e-7	
		2.50000e-1	6.40414e5	3.90372e-7	
2.168	1	1.00000e-2	1624.33362	6.15637e-6	beta-BHC
		5.00000e-2	1.29239e4	3.86880e-6	
		1.00000e-1	3.41054e4	2.93209e-6	
		2.50000e-1	1.94642e5	1.28441e-6	
2.246	1	1.82000e-1	959.13055	1.89755e-4	PCP m.e.
2.326	1	1.00000e-2	4390.62061	2.27758e-6	Lindane
		5.00000e-2	6.35680e4	7.86559e-7	
		1.00000e-1	1.85113e5	5.40211e-7	
		2.50000e-1	6.33512e5	3.94626e-7	
2.368	1	1.00000e-2	3810.67017	2.62421e-6	delta-BHC
		5.00000e-2	5.41755e4	9.22926e-7	
		1.00000e-1	1.79098e5	5.58354e-7	
		2.50000e-1	5.75522e5	4.34388e-7	
2.384	1	2.08000e-1	986.20221	2.10910e-4	2,4,5-T m.e.
2.569	1	5.00000e-2	41.58520	1.20235e-3	Diazinon
2.596	1	1.00000e-1	1616.67700	6.18553e-5	PCNB
2.601	1	7.98000e-1	4760.92773	1.67614e-4	Parathion Methyl
2.639	1	9.99000e-1	1.02025e4	9.79172e-5	Chlorothalonil
2.671	1	5.00000e-2	24.21532	2.06481e-3	Disulfoton
2.722	1	1.89200e-1	28.40828	6.66003e-3	Dinoseb m.e.
2.835	1	5.53500e-1	251.76421	2.19849e-3	Reldan
2.936	1	2.30400e-1	1053.45227	2.18709e-4	Picloram m.e.
2.986	1	5.00000e-2	168.08116	2.97475e-4	Methyl Parathion
3.050	1	1.00200	3852.64917	2.60081e-4	Vinclozolin
3.067	1	1.00000e-2	3353.19849	2.98223e-6	Heptachlor
		5.00000e-2	4.44843e4	1.12399e-6	
		1.00000e-1	1.34210e5	7.45101e-7	
		2.50000e-1	5.42166e5	4.61113e-7	
3.300	1	5.00000e-2	125.14570	3.99534e-4	Malathion
3.319	1	1.00000e-1	1.21988e4	8.19752e-6	Prodiamine
3.364	1	1.56000e-1	3414.49756	4.56875e-5	Chlorpyrifos
3.421	1	1.00000e-2	4408.12842	2.26854e-6	Aldrin

RetTime [min]	Lvl Sig	Amount [ng/ul]	Height	Amt/Height	Ref	Grp Name
		2 5.00000e-2	6.02076e4	8.30460e-7		
		3 1.00000e-1	1.69816e5	5.88873e-7		
		4 2.50000e-1	5.87971e5	4.25191e-7		
3.624	1	5 5.00000e-1	309.82510	1.61381e-3		Dacthal
3.704	1	1 1.00000	311.89117	3.20625e-3		
3.789	1	1 1.00000e-2	3742.91211	2.67172e-6		Fipronil
		2 5.00000e-2	4.59729e4	1.08760e-6		
		3 1.00000e-1	1.33987e5	7.46341e-7		
		4 2.50000e-1	5.15816e5	4.84669e-7		
3.791	1	1 1.00000e-1	496.72495	2.01319e-4		Isofenphos
3.824	1	1 9.88000e-1	2367.17725	4.17375e-4		Captan
3.867	1	1 1.54000	8397.45996	1.83389e-4		Pendimethalin
4.028	1	1 1.96000	261.18826	7.50417e-3		Pirimiphos Methyl
4.063	1	1 1.00000e-2	4995.58008	2.00177e-6		g-chlordane
		2 5.00000e-2	7.46090e4	6.70160e-7		
		3 1.00000e-1	1.95738e5	5.10887e-7		
		4 2.50000e-1	6.40596e5	3.90262e-7		
4.230	1	1 1.00000e-2	4896.07666	2.04245e-6		Endosulfan I
		2 5.00000e-2	6.99753e4	7.14538e-7		
		3 1.00000e-1	1.89318e5	5.28212e-7		
		4 2.50000e-1	6.05361e5	4.12977e-7		
4.260	1	1 1.00000e-2	5996.03906	1.66777e-6		alpha-chlordane
		2 5.00000e-2	9.27057e4	5.39341e-7		
		3 1.00000e-1	2.40732e5	4.15400e-7		
		4 2.50000e-1	6.99314e5	3.57493e-7		
4.448	1	1 1.00000e-2	5544.26660	1.80367e-6		p,p'-DDE
		2 5.00000e-2	9.07221e4	5.51134e-7		
		3 1.00000e-1	2.27711e5	4.39153e-7		
		4 2.50000e-1	7.06875e5	3.53669e-7		
4.481	1	1 1.00000e-2	7179.67969	1.39282e-6		Dieldrin
		2 5.00000e-2	1.15595e5	4.32545e-7		
		3 1.00000e-1	2.90734e5	3.43957e-7		
		4 2.50000e-1	7.75615e5	3.22325e-7		
4.500	1	1 1.00000	0.00000	0.00000		Myclobutanil
4.618	1	1 1.00000e-1	2.62434e4	3.81048e-6		Oxadiazon
4.619	1	1 4.00000e-1	1219.73291	3.27941e-4		Oxyfluorfen
4.649	1	1 1.00000e-2	3986.19043	2.50866e-6		Endrin
		2 5.00000e-2	5.70391e4	8.76592e-7		
		3 1.00000e-1	1.57978e5	6.33000e-7		
		4 2.50000e-1	5.75816e5	4.34166e-7		
4.670	1	1 1.00000e-2	4006.11938	2.49618e-6		Endosulfan II
		2 5.00000e-2	5.87382e4	8.51235e-7		
		3 1.00000e-1	1.70815e5	5.85429e-7		
		4 2.50000e-1	5.87546e5	4.25499e-7		
4.784	1	1 1.00000e-2	4378.14014	2.28407e-6		p,p'-DDD
		2 5.00000e-2	6.55576e4	7.62688e-7		
		3 1.00000e-1	1.86131e5	5.37256e-7		
		4 2.50000e-1	6.29412e5	3.97196e-7		
4.802	1	1 1.00000e-2	3762.21533	2.65801e-6		Endrin Aldehyde
		2 5.00000e-2	5.21170e4	9.59380e-7		
		3 1.00000e-1	1.47239e5	6.79168e-7		
		4 2.50000e-1	5.41579e5	4.61613e-7		
4.966	1	1 1.00000e-2	225.96364	4.42549e-5		Ethion
5.032	1	1 1.00000e-2	764.87964	1.30740e-5		Endosulfan Sulfate
		2 5.00000e-2	7335.86182	6.81583e-6		
		3 1.00000e-1	2.23881e4	4.46666e-6		
		4 2.50000e-1	1.47013e5	1.70053e-6		
5.162	1	1 1.00000e-2	1596.64819	6.26312e-6		p,p' - DDT
		2 5.00000e-2	2.21308e4	2.25929e-6		
		3 1.00000e-1	8.04960e4	1.24230e-6		
		4 2.50000e-1	4.77851e5	5.23176e-7		
5.375	1	1 1.00000e-2	5404.40430	1.85034e-6		Endrin Ketone
		2 5.00000e-2	8.71220e4	5.73908e-7		
		3 1.00000e-1	2.38645e5	4.19032e-7		
		4 2.50000e-1	6.95434e5	3.59488e-7		

RetTime [min]	Lvl Sig	Amount [ng/ul]	Height	Amt/Height	Ref	Grp Name
5.593	1	1.00000e-2	629.91046	1.58753e-5		Methoxychlor
		2 5.00000e-2	5438.36963	9.19393e-6		
		3 1.00000e-1	1.57111e4	6.36493e-6		
		4 2.50000e-1	1.05376e5	2.37246e-6		
5.656	1	1.41500e-1	2847.39185	4.96946e-5		Bifenthrin
5.807	1	5.00000e-2	66.15516	7.55799e-4		Azinphos-methyl
5.937	1	1.00000	0.00000	0.00000		Cyhalothrin I
5.990	1	1.00000	0.00000	0.00000		Cyhalothrin II
6.067	1	9.90000e-1	8787.55176	1.12659e-4		Fenarimol
6.237	1	3.90310e-1	1106.82239	3.52640e-4	1	Permethrin I
6.279	1	6.09690e-1	1728.93262	3.52640e-4	1	Permethrin II
6.433	1	6.97664e-2	1259.81897	5.53781e-5	2	Cyfluthrin I
6.465	1	1.04845e-3	18.93254	5.53781e-5	2	Cyfluthrin II
6.492	1	7.37328e-2	1331.44360	5.53781e-5	2	Cyfluthrin III
6.506	1	9.54524e-2	1723.64868	5.53781e-5	2	Cyfluthrin IV
6.506	1	1.00000	3.02487e5	3.30593e-6		Boscalid
6.535	1	7.86939e-3	20.49156	3.84031e-4	3	Cypermethrin I
6.566	1	4.72971e-2	123.15967	3.84031e-4	3	Cypermethrin II
6.593	1	4.48335e-2	116.74443	3.84031e-4	3	Cypermethrin III&IV
6.714	1	1.00000e-1	1.59947e4	6.25205e-6		DCBP
6.906	1	1.00000e-1	8955.77441	1.11660e-5		Fenvalerate I
6.920	1	1.00000e-1	0.00000	0.00000		Esfenvalerate I
6.971	1	1.00000e-1	8317.79492	1.20224e-5		Fenvalerate II
6.992	1	1.00000e-1	0.00000	0.00000		Esfenvalerate II
7.013	1	1.00000e-1	230.02556	4.34734e-4		Deltamethrin

## More compound-specific settings:

## Compound: Clopyralid m.e.

Curve Type : Linear  
Origin : Forced

## Compound: Dicamba m.e.

Curve Type : Linear  
Origin : Forced

## Compound: 2,4-D

Curve Type : Linear  
Origin : Forced

## Compound: TCMX

Time Window : From 1.774 min To 1.854 min  
Curve Type : Linear  
Origin : Forced

## Compound: alpha-BHC

Time Window : From 2.047 min To 2.087 min  
Curve Type : Power  
Origin : Connected  
Calibration Level Weights:/  
Level 1 : 1  
Level 2 : 1  
Level 3 : 1  
Level 4 : 1

## Compound: beta-BHC

Time Window : From 2.151 min To 2.181 min  
Curve Type : Power  
Origin : Connected

## Compound: PCP m.e.

Curve Type : Linear  
Origin : Forced

## Compound: Lindane

Time Window : From 2.315 min To 2.348 min  
Curve Type : Power  
Origin : Connected

Compound: delta-BHC

Time Window : From 2.353 min To 2.380 min  
Curve Type : Power  
Origin : Connected

Compound: 2,4,5-T m.e.

Time Window : From 2.378 min To 2.396 min  
Curve Type : Linear  
Origin : Forced

Compound: Diazinon

Curve Type : Linear  
Origin : Forced

Compound: PCNB

Time Window : From 2.574 min To 2.625 min  
Curve Type : Linear  
Origin : Forced

Compound: Parathion Methyl

Curve Type : Linear  
Origin : Forced

Compound: Chlorothalonil

Curve Type : Linear  
Origin : Forced

Compound: Disulfoton

Curve Type : Linear  
Origin : Forced

Compound: Dinoseb m.e.

Curve Type : Linear  
Origin : Forced

Compound: Reldan

Curve Type : Linear  
Origin : Forced

Compound: Picloram m.e.

Curve Type : Linear  
Origin : Forced

Compound: Methyl Parathion

Time Window : From 2.971 min To 2.989 min  
Curve Type : Linear  
Origin : Forced

Compound: Vinclozolin

Curve Type : Linear  
Origin : Forced

Compound: Heptachlor

Time Window : From 3.053 min To 3.084 min  
Curve Type : Linear  
Origin : Connected

Compound: Malathion

Curve Type : Linear  
Origin : Forced

Compound: Chlorpyrifos

Time Window : From 3.354 min To 3.382 min  
Curve Type : Linear

Origin : Forced

Compound: Aldrin  
Time Window : From 3.404 min To 3.438 min  
Curve Type : Linear  
Origin : Connected

Compound: Dacthal  
Time Window : From 3.616 min To 3.641 min  
Curve Type : Linear  
Origin : Forced  
Calibration Level Weights:/  
Level 5 : 1

Compound:  
Time Window : From 3.696 min To 3.716 min  
Curve Type : Linear  
Origin : Forced

Compound: Fipronil  
Time Window : From 3.771 min To 3.807 min  
Curve Type : Linear  
Origin : Connected

Compound: Isofenphos  
Time Window : From 3.785 min To 3.802 min  
Curve Type : Linear  
Origin : Forced

Compound: Captan  
Curve Type : Linear  
Origin : Forced

Compound: Pendimethalin  
Curve Type : Linear  
Origin : Forced

Compound: Pirimiphos Methyl  
Time Window : From 4.009 min To 4.048 min  
Curve Type : Linear  
Origin : Forced

Compound: g-chlordane  
Curve Type : Linear  
Origin : Connected  
Calibration Level Weights:/  
Level 1 : 1  
Level 2 : 1  
Level 3 : 1  
Level 4 : 1

Compound: Endosulfan I  
Time Window : From 4.209 min To 4.243 min  
Curve Type : Linear  
Origin : Forced

Compound: alpha-chlordane  
Time Window : From 4.246 min To 4.281 min  
Curve Type : Linear  
Origin : Connected  
Calibration Level Weights:/  
Level 1 : 1  
Level 2 : 1  
Level 3 : 1  
Level 4 : 1

Compound: p,p'-DDE  
Time Window : From 4.426 min To 4.455 min

```

Curve Type           : Linear
Origin               : Connected

Compound: Dieldrin
Time Window          : From 4.465 min To 4.503 min
Curve Type           : Linear
Origin               : Forced

Compound: Oxyfluorfen
Curve Type           : Linear
Origin               : Forced

Compound: Endrin
Time Window          : From 4.633 min To 4.658 min
Curve Type           : Quadratic
Origin               : Connected

Compound: Endosulfan II
Time Window          : From 4.658 min To 4.680 min
Curve Type           : Quadratic
Origin               : Connected

Compound: p,p'-DDD
Time Window          : From 4.775 min To 4.795 min
Curve Type           : Linear
Origin               : Connected

Compound: Endrin Aldehyde
Time Window          : From 4.791 min To 4.815 min
Curve Type           : Quadratic
Origin               : Connected

Compound: Ethion
Time Window          : From 4.941 min To 4.969 min
Curve Type           : Linear
Origin               : Forced

Compound: Endosulfan Sulfate
Time Window          : From 5.023 min To 5.075 min
Curve Type           : Quadratic
Origin               : Connected

Compound: p,p' - DDT
Time Window          : From 5.146 min To 5.191 min
Curve Type           : Linear
Origin               : Connected
Calibration Level Weights:/
Level 1              : 1
Level 2              : 1
Level 3              : 1
Level 4              : 1

Compound: Endrin Ketone
Curve Type           : Quadratic
Origin               : Connected
Calibration Level Weights:/
Level 1              : 1
Level 2              : 1
Level 3              : 1
Level 4              : 1

Compound: Methoxychlor
Curve Type           : Quadratic
Origin               : Connected
Calibration Level Weights:/
Level 1              : 1
Level 2              : 1
Level 3              : 1

```

Level 4 : 1

Compound: Bifenthrin  
Curve Type : Linear  
Origin : Forced

Compound: Azinphos-methyl  
Curve Type : Linear  
Origin : Forced

Compound: Cyhalothrin II  
Curve Type : Linear  
Origin : Forced

Compound: Fenarimol  
Curve Type : Linear  
Origin : Forced

Compound: Permethrin I  
Time Window : From 6.223 min To 6.251 min  
Curve Type : Linear  
Origin : Forced

Compound: Permethrin II  
Time Window : From 6.265 min To 6.291 min  
Curve Type : Linear  
Origin : Forced

Compound: Cyfluthrin I  
Time Window : From 6.423 min To 6.441 min  
Curve Type : Linear  
Origin : Forced

Compound: Cyfluthrin II  
Time Window : From 6.454 min To 6.475 min  
Curve Type : Linear  
Origin : Forced

Compound: Cyfluthrin III  
Time Window : From 6.484 min To 6.500 min  
Curve Type : Linear  
Origin : Forced

Compound: Cyfluthrin IV  
Time Window : From 6.499 min To 6.514 min  
Curve Type : Linear  
Origin : Forced

Compound: Cypermethrin I  
Time Window : From 6.524 min To 6.550 min  
Curve Type : Linear  
Origin : Forced

Compound: Cypermethrin II  
Time Window : From 6.551 min To 6.578 min  
Curve Type : Linear  
Origin : Forced

Compound: Cypermethrin III&IV  
Time Window : From 6.580 min To 6.598 min  
Curve Type : Linear  
Origin : Forced

Compound: DCBP  
Curve Type : Linear  
Origin : Forced

Compound: Deltamethrin

Curve Type : Linear  
Origin : Forced

Group summary :

Group 1 ( Permethrin ) :

Group members:

Permethrin I with retention time 6.237 min  
Permethrin II with retention time 6.279 min

Group Amount Calculation:

Level 1 with amount 1.00000 ng/ul

Group 2 ( Cyfluthrin ) :

Group members:

Cyfluthrin I with retention time 6.433 min  
Cyfluthrin II with retention time 6.465 min  
Cyfluthrin III with retention time 6.492 min  
Cyfluthrin IV with retention time 6.506 min

Group Amount Calculation:

Level 1 with amount 0.24000 ng/ul

Group 3 ( Cypermethrin ) :

Group members:

Cypermethrin I with retention time 6.535 min  
Cypermethrin II with retention time 6.566 min  
Cypermethrin III&IV with retention time 6.593 min

Group Amount Calculation:

Level 1 with amount 0.10000 ng/ul

24 Warnings or Errors (10 first messages follow) :

- Warning : Curve requires more calibration points., (Methomyl)
- Warning : Zero or negative calibration curve slope, (Cyhalothrin II)
- Warning : Curve requires more calibration points., (Boscalid)
- Warning : Overlapping peak time windows at 2.368 min, signal 1
- Warning : Overlapping peak time windows at 2.569 min, signal 1
- Warning : Overlapping peak time windows at 2.596 min, signal 1
- Warning : Overlapping peak time windows at 3.05 min, signal 1
- Warning : Overlapping peak time windows at 3.3 min, signal 1
- Warning : Overlapping peak time windows at 3.789 min, signal 1
- Warning : Overlapping peak time windows at 4.028 min, signal 1

=====  
Peak Sum Table  
=====

\*\*\*No Entries in table\*\*\*  
=====